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MICROSCOPIE A EFFET TUNNEL

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Practical presentation of the STM microscope

The microscope used for this practical work is a commercial STM from Easyscan, sold by the Nanosurf Company.

In general the STM microscope contains three parts: the head, the electronic module and the software control.

- **STM head**

The microscope head is composed by a rigid part where is fixed one or several piezoelectric ceramics. A piezoelectric material has the property to deform itself when it is submitted to an external electrical field. In reverse, a mechanical constraint leads to the apparition of electric charges on its surface. The piezoelectric effect has been discovered by Curie in 1880. In the STM microscope this effect is used to move the tip on the surface with nanometric resolution.

If a voltage is applied between the two faces of a piezoelectric material which are separated by the distance h , an electric field $E = V/h$ appears between the two faces and induces a deformation of the material $\delta h/h = d_{33} \cdot E$ in the electric field direction, where d_{33} represents the piezoelectric constant in the polarisation axis. Typically, the piezoelectric coefficients do not exceed few Å by volt.

The STM microscope used in this practical works (figure 4) is constituted by two elements: the sample holder and the tip holder.


Figure 4 : Easyscan STM head

Tip holder : The tip is fixed to a metallic part (which also allows collecting the tunnelling current), this metallic part is glued on a piezoelectric system, which is constituted by three piezoelectric cubes in order to allow the displacement of the tip in three spatial directions. In presence of a sample and in tunnelling regime, the X and Y movements allow moving the tip across the surface, and the Z direction movement allow keeping the tunnelling current at a constant value.

Sample holder : The sample is glued on a conductive disc. This disc is fixed at the cylinder apex thanks to a magnet placed in the middle of the cylinder. The sample, the metallic disc and the cylinder forms a rigid system. This system can slide on two cylinder tracks in order to approach to the tip. To measure a tunnelling current, the tip-surface distance must be reduce to few nanometers, this approach is realised with an inertial approach motor: The displacement of the cylinder is made by a succession of nano-movements of the support S stuck to a piezoelectric element. During the practical work, we will study the working mode of this approach system. During this approach phase a tunnelling current can be detected, therefore the approach stops immediately, the tunnelling regime is reached.

- ***Electronics module***

- a/ ***Current-voltage amplifier :***

For usual tunnelling voltage (from few mV to few volts), the tunnelling current measured by the STM is relatively low (few pA to few nA). This low value is due to the local character of the measurement. In order to measure with accuracy the I_T current detected by the tip, the current is amplified and converted in a voltage through a current-voltage pre-amplifier located close to the STM head. The figure 5 is a simplified version of this pre-amplifier :


Figure 5 : simplified diagram of the current-voltage pre-amplifier

b/ Control electronics module :

This electronics module which is associated to the PC Pentium 200, is composed by :

- a 16 bits DAC and ADC card : this card converts the analogical signals to digital signal and vice-versa.
- voltage amplifier +/- 12V : these amplifiers are needed to apply voltages at the different piezoelectric elements, from digital ramp generated by the software and converted in analogical signals by the DAC card.

The tunneling current I_T , is converted to a voltage V_S at the output of the preamplifier, and V_S is digitalized by the DAC card. The analogical feedback loop compares V_S to a regulation value fixed by the user, and deduces the error signal ΔZ , which is converted again by the DAC card in analogical signal, which is applied to the piezoelectric element in order to re-find the regulation value. This feedback loop is running continuously in imaging mode, this mode is called constant tunneling current mode. The error signal ΔZ is recorded during the scan X,Y. All the $\Delta Z(X,Y)$ points form the STM image.

This analogical feedback loop is described in figure 6 :


Figure 6 : tunneling current feedback loop

- **Easyscan software**

A detailed description is available during the practical work.

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